GB

Before reading, fold out the illustration page and get to know all of the functions of your unit.

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









GB

Scope of delivery	8
Technical data	8
Correct use	8
Safety notes	8 - 10
Assembly	10 - 11
Training notes	11
Warming up	11 - 12
Use	12
Display function	12 - 13
Stretching	13
Storage, cleaning	13
Disposal	13
Notes on the guarantee and service handling	13 - 14

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
Congratulations!

You have chosen to purchase a high-quality product. Familiarise yourself with the product before using it for the first time.

 **Read the following operating instructions carefully.**

Use the product only as described and only for the given areas of application. Store these operating instructions carefully. When passing the product on to third parties, please also hand over all accompanying documents.

Scope of delivery (figure A)

1 x mobility trainer with resistance screw and display (including battery (1 x 1,5V  LR44))

1 x operating instruction

Technical data

Dimensions: approx. 48 x 30 x 37cm (W x H x D)

Weight: approx. 2.4kg

Material: plastic, steel

Max. user weight: 100kg

Equipment class: HC (for use in private households, precision C) / not suitable for therapeutic training;

Non-speed dependent exercise equipment.

WARNING! This stationary piece of training equipment is not suitable for purposes that require a high level of accuracy.

 Symbol for DC voltage

Date of manufacture (month/year):
02/2020

Correct use

The article is not designed for commercial use. The product is designed for use in the private sphere. This article enables you to train your arms and legs; fitted with a resistance screw in order to increase and maintain mobility. An integrated display shows training advances. Folding legs allow you to store the mobility trainer and assemble it again.

Safety notes

Important: read through these operating instructions carefully and always retain them for future reference.

- The article is not suitable for therapeutic purposes.
- During the exercises sit upright and relaxed with a straight back.
- Exchange defective parts of the article immediately and do not use the article until the repair is complete.
- Please note that incorrect or excessive training may be detrimental to health.
- Place the article on a horizontal, level surface.
- Use an anti-slip underlay on smooth surfaces.
- Observe wear on parts. If the article no longer offers any step resistance then the drive belt is worn and the article needs to be replaced.
- There must be a free area of approx. 0,6m around the training area (figure H, I).


- This device can be used by children aged 14 and above and by persons with reduced physical, sensory, or mental capabilities or with a lack of experience and knowledge only if they are supervised or have received instruction on the safe use of the device and understand the resulting risks. Children may not play with the device. Children may not carry out cleaning or user-maintenance unless supervised.

⚠ Danger to life!

- Never allow children to play with the packaging material unsupervised. There is a risk of suffocation.
- The product is not a piece of climbing or play equipment. Make sure that persons, particularly children, do not stand on the product or pull themselves up on it. The product can tip over.
- The article is not designed to support the entire body weight.
- Children should only use the article under supervision.
- Consult your doctor before beginning use or if you experience any discomfort.

⚠ Risk of injury!

- Do not place your full body weight on the article. It may only be used as described in the instructions.

- Ensure that the article is sufficiently stable before use.
- Place the article on a level surface.
- Be careful of your fingers when unfolding and folding. Fingers may become trapped.
- Ensure that the article is fully extended before use.
- Examine the article for any damage or wear before each use. Replace damages or worn parts on the article immediately as the function and safety could otherwise be impeded.
- The article may only be used for the intended purpose.
- Do not touch the metal parts on the article after use, particularly not the pedal arms, as they heat up during use.

Warning notes - batteries!

- Remove the batteries if they are empty or if the product is not going to be used for an extended period.
- Do not use different battery types or brands, new and used batteries in combination, or batteries with different capacities, as they leak and can thus cause damage.
- Observe the polarity (+/-) when inserting batteries.
- Change all batteries at the same time and dispose of the old batteries properly.
- Warning! Batteries may not be charged or reactivated in any way. They may not be dismantled, thrown into fire, or short-circuited.
- Always keep batteries out of the reach of children.
- Do not use rechargeable batteries!
- The cleaning and user-maintenance may not be carried out by children without supervision.
- Clean the batteries and the device contacts as necessary and before inserting the batteries.

- Do not expose the batteries to extreme conditions (e.g. radiators or direct sunlight). Otherwise there is an increased risk of leaks.
- Batteries can endanger life if swallowed. Keep batteries out of the reach of small children for this reason. Seek medical help immediately if a battery is swallowed.

Danger!

- Handle a damaged or leaking battery with extreme care and dispose of it properly and immediately. Wear gloves to handle the battery.
- If you come into contact with battery acid then wash the affected area with soap and water. If battery acid gets in your eye rinse it out with water and seek medical help immediately!
- Do not short-circuit the terminals.

Assembly

Assembling the article

1. Remove the packaging material.
2. Unfold both legs until you hear an audible click from the safety pins (figure B).
3. Put the article up in front of you and turn the resistance screw clockwise to increase resistance and counter-clockwise to decrease resistance (figure C).

Note: '+' and '-' are marked on the resistance screw to adjust the resistance. '+' stands for increase and '-' stands for reduction of the resistance.

Folding the article (figure D)

1. Turn the article onto its side.
2. Pull on the locking ring on each leg to release the locking pin, then fold the legs in one at a time.

Adjusting the pedal strap (figure E)

The pedal steps are already assembled on the pedals. The strap can be set to one of three widths.

Note: You will require a corresponding screwdriver.

1. Unscrew the two screws with the correct screwdriver and remove them from the pedal and the end of the strap.
2. Select the holes on the end of the strap for the desired width and insert the two screws in the holes on the end of the strap and the pedals.
3. Tighten the screws with the correct screwdriver. Set the widths on both pedals.

Changing the battery

CAUTION! Follow these instructions in order to avoid mechanical and electrical damage.

Danger!

One LR44 battery is included in the article. If the battery is no longer working then you can change it.

1. Lift the display on the side and remove it carefully (figure F). The battery is on the back of the display.
2. Carefully slide the spent battery out of the battery holder and exchange it for a new battery (figure G).

Note: ensure that the battery is completely inside the mounting.

3. Insert the cable connections and the display back into the casing.

Note: please ensure that the display is replaced correctly into the casing.

Training notes

The following exercises only represent a selection. You can find further exercises in the relevant literature.

- Warm up well before each training session.
- Do not train if you feel ill or unwell.

- When starting training, get instruction on the correct execution of the exercises from an experienced physiotherapist if possible.
- Only exercise as long as you feel well and/or you are able to carry out the exercises correctly.
- Wear comfortable sports clothing and trainers.
- Determine the frequency and intensity of the exercises yourself. Start slowly 2 to 3 times a week for 10 minutes at a time, and increase the frequency and the exercise intensity gradually. The more frequently and regularly you carry out the exercises, the fitter and better you will feel.

Caution! Avoid training too intensely!

- When you first begin training, 2-3 minutes are sufficient per exercise. If you train every day the duration can be increased to 5-10 minutes after approximately one week. The maximum training time should, however, not exceed 1 hour.
- Take sufficient breaks between exercises and drink sufficient fluids.

Caution!

Cease exercising immediately if you experience discomfort or feel unwell and consult your doctor.

Warming up

Take sufficient time before each training session to warm up. The following includes descriptions of some simple exercises for this purpose. Repeat these exercises 2 to 3 times each.

Warming up neck muscles

1. Turn your head slowly to the left and to the right.
2. Repeat this movement 4-5 times.
3. Circle your head slowly, first in one direction and then in the other.

Warming up arms and shoulders

1. Circle both your shoulders forwards at the same time.
2. Change direction after one minute.
3. Then pull your shoulders up to your ears and let your shoulders fall again.

4. Circle your right and left arms alternately forwards and then after one minute backwards.

Important: Don't forget to continue breathing steadily throughout.

Warming up leg muscles

1. Stand on one leg and lift the other leg with your knee bent until it is approximately 20cm from the floor.
2. First circle the lifted foot in one direction and then change direction after several seconds.
3. Then change legs and repeat the exercise.
4. Lift your legs in turn and take several steps in place. Make sure that you lift your legs only as far as you are able to maintain good balance.

Use

Place the training article in front of you. Start with easy training and low resistance. Release the resistance screw to achieve this. The tighter you turn the resistance screw, the more difficult it is to turn the pedals. Begin your training in the warm-up phase with slow movements and low resistance. Increase the speed and the resistance after 5-10 minutes. Reduce the speed and the resistance again at the end of your training session. Carry out some stretching exercises for the corresponding muscles after your training session.

Leg training (figure H)

For this exercise the chair must be placed in front of the article so that your knees always remain slightly bent in the extended phase. Ensure that you sit on the chair without support and do not lean back.

1. Place the article on a level, stable and non-slip area of floor, and sit upright and with a straight back on a chair in a relaxed position in front of the article.

Note: ensure that the chair is stable and comfortable.

2. Place your feet on the pedals so that pressure is placed on the centre of the soles.
3. Tip your pelvis forwards and tense your abdominal muscles. Keep your abdominal muscles tensed throughout the exercise.

4. Keep your upper body upright and pull your shoulder blades to your spine. Your head should form the extension of your spine.
5. 'Cycle' with your feet and avoid evasive movements inwards and outwards.

Arm training (figure I)

For this exercise the article must stand in a position so that the movement takes place at the height of your breastbone or under your shoulders. Ensure that your elbows are slightly bent in the extended phase.

1. Place the article in front of you on a level, stable and non-slip surface (e.g. at table height), and sit upright with your back straight in a relaxed position on a chair.
2. Place your feet on the floor shoulder-width apart and tip your pelvis forwards.
3. Tense your abdominal muscles, straighten your upper body and pull your shoulder blades towards your spine. Keep your abdominal muscles tensed throughout the exercise. Keep your shoulders down.
4. Place the palms of your hands on the pedals so that the pedal strap is above your hand.
5. 'Cycle' with your hands. Your head should form the extension of your spine. Ensure that you do not overextend your head.

Display function

The display has various modes that you can select by pressing the button:

- SCAN: This mode displays the various display modes alternately.
- TIME: This mode displays the duration of your current training session.
- CNT: This mode displays the rounds completed in the current training session.
- CAL: This mode displays the calories used during your current training session.

Note: calorie consumption is only an approximate reference value.

- R.P.M.: This mode displays the speed in rounds per minute.

Note: the display switches off automatically after approx. four minutes when not in use.

Set value to zero

Press on the button for approximately four seconds to reset the value to 0.

Stretching

Take sufficient time after each training session to stretch. The following illustrates a simple exercise for this purpose. You should carry out the stretching exercise 1 times on each side for 10-15 seconds.

Neck muscles

1. Stand in a relaxed position. Pull your head gently to the left and to the right using one hand. This exercise stretches the sides of your neck.

Arms and shoulders

1. Stand up straight with your knees slightly bent.
2. Put your right arm behind your head until your right hand is between your shoulder blades.
3. Hold your right elbow with your left hand and pull it backwards.
4. Change sides and repeat the exercise.

Leg muscles

1. Stand up straight and lift one foot off the ground.
2. Slowly circle your foot, first in one direction and then in the other.
3. Change feet after a while.

Important: Please ensure that your thighs stay parallel next to one another. Push your hips forward and keep your upper body upright.

Storage, cleaning

Fold the legs away as described in the section 'folding the article'. Store the article in a dry place.

Parts that are susceptible to wear: if the article is no longer offering any step resistance then the drive belt is worn. The article must be replaced. Clean the article with a damp cloth and a mild cleaning agent.

Caution! Do not allow moisture to penetrate the display.

Disposal

In the interest of protecting the environment, do not throw your product into the household waste once you are finished with it, but rather take it to a speciality disposal facility. Find out about collection sites and their hours of operation from your local authority.

Defective or used batteries must be recycled pursuant to Directive 2006/66/EC and changes to it. Return batteries and/or the product via the collecting sites provided. Packaging materials, such as plastic bags, must be kept away from children. Keep the packaging materials out of reach of children.

Environmental damage caused by improper disposal of batteries! Batteries may not be disposed of with household waste. They may contain toxic heavy metals and are subject to hazardous waste treatment. The chemical symbols of the heavy metals are as follows: Cd = cadmium, Hg = mercury, Pb = lead. Dispose of used batteries at a municipal collection site.

Dispose of the product and packaging in an environmentally friendly manner.

Note the label on the packaging materials when separating waste, as these are labelled with abbreviations (a) and numbers (b) with the following meaning: 1-7: plastics/20-22: paper and cardboard/80-98: composite materials.

The product and the packaging materials can be recycled, dispose of them separately for better treatment of waste. The Triman logo only applies to France. Find out how to dispose of the used product from your municipal or city administration.

Notes on the guarantee and service handling

The product was produced with great care and under constant supervision. You receive a three-year warranty for this product from the date of purchase. Please retain your receipt.

The warranty applies only to material and workmanship and does not apply to misuse or improper handling. Your statutory rights, especially the warranty rights, are not affected by this warranty.

With regard to complaints, please contact the following service hotline or contact us by e-mail.

Our service employees will advise as to the subsequent procedure as quickly as possible.

We will be personally available to discuss the situation with you.

Any repairs under the warranty, statutory guarantees or through goodwill do not extend the warranty period.

This also applies to replaced and repaired parts.

Repairs after the warranty are subject to a charge.

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